Docket No: 1743/227

IN THE CLAIMS:

1. Canceled.

2. (New) A specimen observing method, comprising the steps of:

irradiating a specimen supplied with a negative voltage with an electron beam to

generate secondary electrons from the specimen,

deflecting and detecting the generated secondary electrons to obtain a specimen

image on the basis of the detected secondary electrons, and

adjusting the negative voltage so as to cancel charge-up of the specimen.

3. (New) A specimen observing method according to claim 2, wherein the step of adjusting the negative voltage is so performed that the specimen image has a given image

contrast.

4. (New) A specimen observing method according to claim 2, wherein the step of

adjusting the voltage is performed while monitoring an output of a secondary electron detector

detecting the generated secondary electrons so that the output shows a maximum.

5. (New) A specimen observing method, comprising the steps of:

irradiating a specimen supplied with a negative voltage with an electron beam to

generate secondary electrons from the specimen,

deflecting and detecting the generated secondary electrons to obtain a specimen

image on the basis of the detected secondary electrons, and

adjusting the negative voltage while monitoring an output of a secondary electron

detector detecting the secondary electrons to determine a value of the negative voltage at which

the output has a maximum value.

2

PATENT

Serial No: 10/699,793

Docket No: 1743/227

6. (New) A scanning electron microscope which comprises:

an electron source emitting a primary electron beam,

a focusing lens focusing the primary electron beam,

a deflector deflecting the primary electron beam toward the specimen beam to irradiate a specimen with the focused primary electron beam so as to generate secondary

electrons from the specimen;

a detector detecting the secondary electrons;

a negative voltage supplying source supplying the specimen with a negative

voltage, and

a controller adjusting the negative voltage to determine a value thereof at which

the output shows a maximum.

7. (New) A scanning electron microscope according to claim 6, wherein the secondary electron detector comprises a secondary electron multiplier and the controller adjusts the negative voltage while monitoring an output of the secondary electron multiplier.